


Probe Software

Probe Software, Inc.
885 Crest Drive
Eugene, OR 97405 USA
(541) 343-3400
sales@probesoftware.com
www.probesoftware.com




Software for MicroAnalysis

Probe for EPMA

Probe Image

PictureSnapApp

TDI (time dependent intensity) Scanning for
Quantitative X-ray Mapping of Beam Sensitive Samples



Using fast replicate scans we are able to extrapolate x-ray map intensities to zero time, on a pixel by pixel basis, for quantitative correction of alkali ion migration (or carbon contamination) and other beam damage artifacts.